





## One size, fits left and right



Telescopic-bar extension system to enable height adjustment from 46 to 69 cm.

## Non-slip-textured breathable pads

Padded non-slip-textured medial-lateral plates that prevent the orthosis from slipping or rotating around the leg.

#### Joint with single locking

Locking with single button. Easy and comfortable: 0°, 15° and 30°



Introducing the new optimiced e-sense telescopic ROM knee brace

for controlling and protecting the flexion/extension locking and movement of the knee in pre/post-surgical or post-injury situations.

## maximum Control, optimum recovery



# Wide range of flexion/extension control

Flexion control from 0° to 120° in increments of 10° and extension control from 0° to 90° also in increments of 10°

## Tool-less adjustment

Simple assembly, no tools required for adjustment.

#### Secure fastening

Straps to ensure secure fastening and perfect adjustment to the degree of inflammation of the knee.



## maximum **C**ontrol, optimum recovery





Controlled and protected movement



**Telescopic ROM knee brace** *RS-3000* 

P. Industrial №1 Calle F, 15 28938 Móstoles - Madrid WWW. Drim. ⊖S

